

Digital Dry Bath Incubator (Interchangeable Blocks)

MODEL LT-442



INTRODUCTION

Digital Dry Bath Incubator (Interchangeable Blocks) is a unique instrument for preservation & reaction of samples, Amplification of DNA, pre-denaturation of the electrophoresis, Blood serum coagulation in molecular biology, histology, clinical, environmental and industrial settings and incubation at regulated Temperature of the Dry Bath.

TECHNICAL SPECIFICATION :

Display	16 X 2 Line Alpha Numeric LCD
Display Format	Menu driven
Hole (Interchangeable Blocks)	25 Holes for 12.5 mm tubes
Setting	By Soft Touch Keys
Heating Devices	Test Tube Dry Bath
Temp. Sensing Devices	RTD/ PT100/ Digital Semiconductor
Range	5°C above ambient to 150°C
Resolution	1°C
Accuracy	± 1°C
Calibration	Manual (if required)
Block size	24 x 12 mm Tube holders.
Block Material	Rock solid Aluminum block.
Power Supply	230 V ±10% AC 50 HZ

LABTRONICS®
Always a step ahead in technology...

CE
ISO9001

एन एस आई सी
NSIC
CERTIFIED