

LT-10/ LT-11

DIGITAL pH METER

PRODUCT FEATURES

- ✓ Compact design
- ✓ Bright LED display
- ✓ Accurate and stable result
- ✓ 3 point calibration (LT-10)
- ✓ 2 point calibration (LT-11)

ACCESSORIES

- ✓ Combination pH Electrode
- ✓ Buffer Tablets 4 pH & 7 pH
- ✓ Operation Manual
- ✓ Buffer Bottles (4 pH & 7 pH)
- ✓ Dust Cover
- ✓ Temperature Probe (in Model LT-10)
- ✓ Rotational pH Electrode Stand



TECHNICAL SPECIFICATIONS

Parameter	pH	mV
Range	0.00 to 14.00 pH	0 to 1999.9 mV
Resolution	0.01pH	1 Mv
Accuracy	±0.01 pH	± 1mV
Display	3½ Digit LED	
EMP Compensation	Manual : 0 to 100 1C (ATC in Model LT - 10)	
Input Impedance	> 10 ohms	
Slope Control	80 to 120%	
Recorder Output	0 to 10 mV/pH Adjustable 0 to 10 mV/100 mV Adjustable	
Power	230 V ±10% AC, 50Hz	
Dimensions	275 x175x 76 mm (LxBxH)	
Sensor	Combined pH Electrode	
Weight	2kg	