

LT-101

DIGITAL pH METER

INBUILT STIRRER

PRODUCT FEATURES

- ✓ Inbuilt stirrer
- ✓ Bright LED display
- ✓ Stable and accurate result
- ✓ Auto temperature compensation

ACCESSORIES

- ✓ Combination pH Electrode
- ✓ Buffer Tablets 4 pH & 7 pH
- ✓ Operation Manual
- ✓ Buffer Bottles (4 pH & 7 pH)
- ✓ Dust Cover
- ✓ Temperature Probe
- ✓ Rotational pH Electrode Stand
- ✓ Magnetic Capsule



TECHNICAL SPECIFICATIONS

Parameter	pH	mv
Range	0.00 to 14.00 pH	0 to 1999.9 mV
Resolution	0.01pH	1 Mv
Accuracy	± 0.01 pH	± 1mV
Display	3 1/2 Digit LED	
Temp Compensation	Manual : 0 to 100 1C (ATC in Model LT - 10)	
Input Impedance	> 10 ohms	
Slope Control	80 to 120%	
Recorder Output	0 to 10 mV/pH Adjustable 0 to 10 mV/100 mV Adjustable	
Power	230 V ±10% AC, 50Hz	
Dimensions	275 x175x 76 ctn (LxBxH)	
Sensor	Combined pH Electrode	
Stirrer	inbuilt Stirrer (with speed variable)	