

# LT-109 MICROPROCESSOR MELTING POINT APPARATUS

AUTOMATIC WITHOUT SILICON OIL

## PRODUCT FEATURES

- ✓ 128x64 Graphical LCD display with backlight
- ✓ Icon-based navigation menu for ease of operation
- ✓ Selectable heating modes (slow, medium and Rapid)
- ✓ Storage up to 1000 results
- ✓ Quick Cooling with Integrated fans
- ✓ Audio and Visual Indicators for System Status
- ✓ Quick Cooling with integrated Fans

## ACCESSORIES

- ✓ Capillaries 1 Box
- ✓ Cleaning tool 1 Nos.
- ✓ USB Cable 1 Nos.
- ✓ Software 1 Nos.



## TECHNICAL SPECIFICATIONS

No of sample	1 sample at a time
Temperature range	Ambient to 300 °C
Temperature accuracy	±1 °C at 200 °C and ±2.5 °C at 275 °C
Temperature sensor	Platinum resistance PT-100
Selectable heating mode	Slow, medium and rapid
Display	128x64 graphical LCD display
No. of key	4
Memory	Storage upto 1000 results
Power	230V ±10% AC 50hz
Dimension	300x200x110mm (approx)
Weight	3kg



**LABTRONICS®**  
Always a step ahead in technology...

All Designs And Specifications are subject to change without prior notice.