

LT-110 DIGITAL DRY MELTING POINT APPARATUS


VISUAL - WITHOUT SILICON OIL

PRODUCT FEATURES

- ✓ Furnace based dry heating
- ✓ Highly accurate and stable
- ✓ Excellent view angle
- ✓ Fast heating

ACCESSORIES

- ✓ Operation Manual
- ✓ Capillary Tubes
- ✓ Power Cable
- ✓ Test Tube for Boiling Samples



TECHNICAL SPECIFICATIONS

Display	3 Digit LED
Display Format	Menu driven
Setting	By Soft Touch Keys
Heating Devices	Dry Furnace
Temp. Sensing Devices	RTD/ PT100
Range	5°C above ambient to 350°C
Resolution	1°C
Accuracy	Depend on manual viewing/ operation
Calibration	Manual
Power Supply	230 V \pm 10% AC 50 HZ
Dimension	355x340x450mm (approx)
Weight	6kg