

# LT-113P PORTABLE HAEMOGLOBIN METER

## PRODUCT FEATURES

- ✓ Compact body design
- ✓ Accurate and stable result
- ✓ 1000 test results

## ACCESSORIES

- ✓ HB meter
- ✓ Canister of test strips
- ✓ Code chip
- ✓ Control chip
- ✓ Capillary transfer tube/dropper
- ✓ Lancing device
- ✓ Sterile lancets
- ✓ AAA batteries



## TECHNICAL SPECIFICATIONS

Testing Strips	10 strips
Measuring Range	5-25 gm/dl
Memory	1000 test results
Calibration	Auto calibration