

LT-114 DIGITAL PHOTO COLORIMETER

AUTO ZERO

PRODUCT FEATURES

- ✓ Auto zero facilities
- ✓ Accurate and stable result
- ✓ Compact design
- ✓ Bright LED display
- ✓ Long life silicon photodiode

ACCESSORIES

- ✓ Glass test tubes 5 Nos.
- ✓ Instruments cover 1 Nos.
- ✓ Operating Manual 1 Nos.
- ✓ Test tube stand 1 Nos
- ✓ Dust cover 1 Nos



TECHNICAL SPECIFICATIONS

Filters	8 filters
Accuracy	±0.02 O.D.
Reproducibility	±0.01 O.D.
Stability	±0.02 O.D per hour
Light Source	3.5V LED Bulb
Key Board	1 Soft Touch Key
Amplifier	I.C. Based
Detector	Selenium photocell/Photodiode
Warm-up Time	5 Minutes.
Sample System Tubes	12 mm dia matched glass test
Power Source	230V ±10 % 50Hz.AC.
Dimension	250x200x100mm approx
Weight	3kg approx