

LT-12 DIGITAL PHOTO COLORIMETER

PRODUCT FEATURES

- ✓ Bright LED display
- ✓ Accurate and stable result
- ✓ Set zero with potentiometer adjustment
- ✓ Long life silicon photodiode

ACCESSORIES

- ✓ Glass test tubes 5 Nos.
- ✓ Instruments cover 1 Nos.
- ✓ Operating Manual 1 Nos.
- ✓ Test tube stand 1 Nos
- ✓ Dust cover 1 Nos



TECHNICAL SPECIFICATIONS

Filters	8 Filters
Accuracy	± 0.02 O.D.
Reproducibility	± 0.01 O.D.
Stability	± 0.02 O.D. Per Hour
Light Source	3.5V LED Bulb
Key Board	Potentiometer Adjustment
Amplifier	I.C. Based
Detector	Selenium Photocell/photodiode
Warm-up Time	5 Minutes.
Sample System Tubes	12mm Dia Matched Glass Test
Power Source	230V \pm 10% 50hz.ac
Dimension	280x260x160mm approx
Weight	3kg