

LT-1261 MICROPLATE ELISA WASHER

PRODUCT FEATURES

- ✓ Large LCD display, intelligent software, easy to operate
- ✓ Strong washing function, bottom washing and two positions aspiration make completed washing effect
- ✓ 8-channel or 12 channel optional
- ✓ Needle positions adjustable and up to 10 different needle positions pre-settable
- ✓ Pressure, volume & time adjustable, washing times, soap time pre-settable
- ✓ Pressure, volume & time adjustable, washing times, soap time pre-settable, whole plate or single strip reading
- ✓ Self calibration and trouble alarm function
- ✓ Large memory to store up to 100 program



ACCESSORIES

- ✓ 4 wash bottle(2500ml)
- ✓ Connecting tube set
- ✓ Cleaning wire
- ✓ Main lead
- ✓ Instruction manual



LCD
DISPLAY



ACCURATE
& STABLE



SOFTWARE
CONTROL SYS.



EASY
OPERATION



RS232
CONNECTIVITY

TECHNICAL SPECIFICATIONS

Plate	48/96 well and trips (U.V. Flat bottom)
Residual Volume	≤ 2µl
Dispensing volume	Adjustable from 50 to 900 µl
Wash cycles	1-16 optional
Mode of wash	From 1 to 12, row and rank, selectable
Soak time	Up to 990 seconds adjustable
Shake time	Up to 990 seconds adjustable
Display	Large LCD
Channel	3 channels for washing liquid, 1 channel for waste liquid
Memory	up to 100 programs can be stored
Interface	RS232
Dimension	450(L)mm x 310(W)mm x 170 (H)mm
Net Weight	7.5kg
Power	AC220V ± 10%, 50/60Hz; AC110V ± 10%, 50/60Hz