

LT-15 DIGITAL TDS METER

PRODUCT FEATURES

- ✓ Bright LED display
- ✓ Stable and accurate result
- ✓ Compact and accurate
- ✓ Easy to operate

ACCESSORIES

- ✓ TDS cell 1.0
- ✓ Main lead
- ✓ Rotational stand
- ✓ Instruction manual



TECHNICAL SPECIFICATIONS

Display	3½ Digit LED
Measurement	TDS
Cond. ranges	Five Ranges
Range Selection	Automatic in 5 Ranges
Temp. Range Accuracy	±5 to 110°C
Cond.	±0.5% FS ± 1DIGIT
Temperature Compensation	Auto: 0 to 50 °C Manual: 0 to 50 °C
Cell Constant	Adjustable on Digital Display
Measuring Cell	Platinum DIP Type
Cell Const Measurement	Auto
Power	230V ± 10% AC
Dimensions	270x202x80 mm (lxbxh)
Weight	0.8 Kg approx