

# LT-151 DIGITAL TDS METER

## INBUILT STIRRER



### PRODUCT FEATURES

- ✓ Inbuilt stirrer
- ✓ Bright LED display
- ✓ Stable and accurate result
- ✓ Easy to operate

### ACCESSORIES

- ✓ TDS cell 1.0
- ✓ Rotational TDS cell stand
- ✓ Magnetic capsule
- ✓ Instruction manual
- ✓ Dust cover



### TECHNICAL SPECIFICATIONS

Display	3½ Digit LED
Measurement	TDS
Cond. Ranges	Five Ranges
Range Selection	Automatic 5 Ranges
Temp. Range	±5 to 110°C
<b>ACCURACY</b>	
Cond.	±0.5% FS ± 1DIGIT
Temp	±0.1°C ± 1DIGIT
Stirrer	Inbuilt Stirrer (With Speed Variable)
Temperature Compensation	Auto: 0 to 50 °C Manual: 0 to 50 °C
Cell Constant	Adjustable on Digital Display
Measuring Cell	Platinum DIP Type
Cell Const	Auto
<b>MEASUREMENT</b>	
Power	230V ± 10% AC
Dimensions	270x202x80 mm (LxBxH)
Weight	0.8 kg approx
Accessories	TDS Cell, Magnetic Capsule, Operation Manual, Dust Cover, Cell Stand,