

LT-26/ LT-16 DIGITAL CONDUCTIVITY METER

PRODUCT FEATURES

- ✓ Bright LED display
- ✓ Automatic temperature compensation (Model LT-26)
- ✓ Stable & accurate result
- ✓ Compact design

ACCESSORIES

- ✓ Conductivity cell 1.0
- ✓ Temperature probe (LT-26)
- ✓ Rotational Conductivity cell stand
- ✓ Instruction manual
- ✓ Dust cover



TECHNICAL SPECIFICATIONS

Display	3½ Digit LED
Measurement	Conductivity
Cond. Ranges	Five Ranges
Range Selection	Automatic in 5 Ranges
Accuracy	
Cond.	± 0.5% FS ± 1DIGIT
Temperature Compensation	Auto 0 to 50°C (in model LT-26 only) manual 0 to 50°C in (Model LT-16)
Cell Constant	Adjustable on Digital Display
Measuring Cell	Platinum DIP Type
Cell Const. Measurement	Auto
Power	230V ± 10% AC
Dimensions	270x202x80 mm (LXBXH)
Weight	0.8 Kg approx.
Temperature probe	For model LT-26