

LT-160 DIGITAL BOMB CALORIMETER

PRODUCT FEATURES

- ✓ Manual measurement & calculation.
- ✓ Initial temperature rise & fall. display
- ✓ Accurate result



ACCESSORIES

- ✓ Steel drum
- ✓ Bomb with sample holder
- ✓ Stirrer
- ✓ Temperature Probe
- ✓ Nichrome Wire
- ✓ Pallet press
- ✓ Cotton thread
- ✓ Safely valve for high pressure
- ✓ Oxygen pressure gauge
- ✓ Benzoic acid powder
- ✓ Tool kit
- ✓ Operation manual
- ✓ Power supply cable (Main lead)



TECHNICAL SPECIFICATIONS

| | |
|------------------------|-----------------|
| Display | 3 Digit LED |
| Temperature accuracy | 0.1°C |
| Temperature Resolution | 0.1°C |
| Temperature Detector | RTD |
| Temp. Range | RT. to 99.9°C |
| Calculate | Manual |
| Power supply | 230V ± 10%/50Hz |