

# LT-161 DIGITAL CONDUCTIVITY METERS

## INBUILT STIRRER

### PRODUCT FEATURES

- ✓ Inbuilt stirrer
- ✓ Bright LED display
- ✓ Stable and accurate result
- ✓ Easy to operate

### ACCESSORIES

- ✓ Magnetic capsule
- ✓ Conductivity cell 1.0
- ✓ Rotational conductivity cell stand
- ✓ Instruction manual
- ✓ Dust cover



### TECHNICAL SPECIFICATIONS

Display	3½ Digit LED
Measurement	Conductivity
Cond. Ranges	Five Ranges
Range Selection	Automatic 5 Ranges
Accuracy	
Cond.	±0.5% FS ± 1DIGIT
Temperature Compensation	Manual 0 to 50°C
Stirrer	Inbuilt Stirrer (With Speed Variable)
Cell Constant	Adjustable on Digital Display
Measuring Cell	Platinum DIP Type
Cell Const. Measurement	Auto
Power	230V ± 10% AC
Dimensions	270x202x80 mm (LXBXH)
Weight	0.8 Kg approx.
Accessories	Conductivity Cell, Magnetic Capsule, Operation Manual, Dust Cover, Cell Stand,