

# LT-175 MICROPROCESSOR AUTO LEAK TEST APPARATUS

## 12" DESICCATOR

### PRODUCT FEATURES

- ✓ Programmable Time With Digital Display
- ✓ Auto Stop After Completion Of Time With Alarm
- ✓ Precise Control Of Vacuum Pump
- ✓ Heavy Duty Oil-Free Vacuum Pump
- ✓ 1 mm Thick SS-304 Cabinet

### ACCESSORIES

- ✓ Instruction manual
- ✓ 12 Inch vacuum desiccator
- ✓ Connecting pipes



LED DISPLAY



ACCURATE & STABLE



SS 304 CABINET



AIR CIRCULATION



EASY OPERATION

### TECHNICAL SPECIFICATIONS

Type	GMP Model, In-built Test Unit
Time Control & Range	Digital, Programmable, 00.01 sec to 99min59sec
Test Vacuum Range	Adjustable up-to 650 mmHg
Vacuum Control Range	Digital Display
Vacuum control Release	By Single Front Panel Knob
Audio Alarm at Completion	Yes
Vacuum Pump	Oil free, MV 27" ± 5%, 1/4HP Motor
Test Unit Size	200mm, 250mm & 300mm (any one) (Specify in Order)
Power	230V ± 10% AC 50hz
Dimension	600x450x200mm (approx)
Weight	8kg (approx)



**LABTRONICS®**  
Always a step ahead in technology...

All Designs And Specifications are subject to change without prior notice.