

LT-181 DISSOLVED OXYGEN METER

INBUILT STIRRER

PRODUCT FEATURES

- ✓ Inbuilt stirrer
- ✓ Bright LED display
- ✓ Compact design
- ✓ Accurate and stable result

ACCESSORIES

- ✓ Gold/silver D.O Probe
- ✓ D.O membrane
- ✓ Rotational D.O probe stand
- ✓ Magnetic capsule
- ✓ Dust cover
- ✓ Instruction manual



TECHNICAL SPECIFICATIONS

Range	0 to 20.0 ppm
Resolution	0.1 ppm
Accuracy	0.1ppm \pm 1 count
Temp. Compensation	0°C to 50°C manual
Compensation	\pm 3%
D.O. Sensor	amperometric gold/silver membrane type.
Stirrer	Inbuilt stirrer with speed control
Display	3½ digit LED display
Power	230V A.C 50 Hz
Weight	2 Kg



LABTRONICS®
Always a step ahead in technology...

All Designs And Specifications are subject to change without prior notice.