

LT-21 DIFFERENTIAL CELL COUNTER




PRODUCT FEATURES

- ✓ Facility to count 8 type of cells
- ✓ 16x2 line alphanumeric LCD display
- ✓ Automatic percentage (%) calculation
- ✓ Individual cell counting facility
- ✓ Audible beep facility

ACCESSORIES

- ✓ Dust cover
- ✓ Instruction manual



TECHNICAL SPECIFICATIONS

Cell Counting Key	16 keys
Total Keys	16 keys (8 cell keys and 3 function keys)
Body	Metal body
Display	16x2 line LCD display
Mains Input	230±10% ac, 50 hz
Weight	1.5 kg (approx)
Dimensions	250x230x150 mm (LXBXH)