

# LT-22 DIGITAL POTENTIO- METER

INBUILT STIRRER

## PRODUCT FEATURES

- ✓ Inbuilt stirrer
- ✓ 3½Digit LED display
- ✓ Highly accurate and stable
- ✓ Compact design

## ACCESSORIES

- ✓ Calomel electrode
- ✓ Glass electrode
- ✓ Platinum electrode (optional)
- ✓ Silver electrode (optional)
- ✓ Magnetic capsule
- ✓ Rotational electrode stand



LED DISPLAY



ACCURATE & STABLE



INBUILT STIRRER



STABLE DESIGN



EASY OPERATION

## TECHNICAL SPECIFICATIONS

Display	3½ digit bright red LED display
Range	0 to $\pm 1999$ mV
Resolution	1 mV
Accuracy	$\pm 1$ mV $\pm 1$ digit
Operating Temperature	10°C to 50°C
Inbuilt Stirrer	Available
Power Supply	230V $\pm 10\%$ AC, 50 Hz
Dimensions	270X 248 X 101 mm (LxBxH)
Weight	2 kg.