

LT-23 DIGITAL pH-COND. TEMPERATURE METER

PRODUCT FEATURES

- ✓ 3 parameters in single instrument (pH/Conductivity/temperature)
- ✓ Compact design
- ✓ Bright LED display
- ✓ Accurate and Stable

ACCESSORIES

- ✓ Conductivity cell 1.0
- ✓ Combination pH Electrode
- ✓ Buffer Tablets 4 pH & 7 pH
- ✓ Operation Manual
- ✓ Dust Cover
- ✓ Rotational pH Electrode Stand
- ✓ Temperature probe



TECHNICAL SPECIFICATIONS

pH MEASUREMENTS

pH	0 to 14 pH
Resolution	0.01 pH
Accuracy	$\pm 0.01 \text{ pH} \pm 1 \text{ digit}$

CONDUCTIVITY METER

Range	0 to 2 & 20mMhos/cm
Accuracy	$0.5\% \pm 1 \text{ digit}$

TEMPERATURE MEASUREMENT

Range	0° to 100°C
Resolution	0.1 ° C
Accuracy	0.2 % + 1 digit
Display	3½ digit LED display
Power	230V $\pm 10\%$ AC 50hz
Dimensation	275x125x175mm approx
Weight	3kg