

# LT-25 DIGITAL TDS- CONDUCTIVITY METER

## PRODUCT FEATURES

- ✓ 2 parameters in single instrument (TDS/Conductivity)
- ✓ Bright LED display
- ✓ Accurate and stable
- ✓ Easy to operate

## ACCESSORIES

- ✓ TDS/COND. Cell ( Dip Type)
- ✓ Rotational stand
- ✓ Main lead
- ✓ Operation Manual
- ✓ Dust Cover



LED DISPLAY



ACCURATE & STABLE



STABLE DESIGN



24 HOURS WORKING



EASY OPERATION

## TECHNICAL SPECIFICATIONS

|                          |   |
|--------------------------|---|
| Range                    | COND. 200 $\mu$ /s to 1000 ms in 5 ranges.                          |
| TDS                      | - 200ppm to 1000 ppt in 5 ranges.<br>All ranges start from zero "0" |
| Accuracy                 | $\pm 0.3\%$ F.S. 200 Micro Mhos/cm to 200m/                         |
| Mhos/cm.                 | $\pm 1\%$ FS. In last range 1000m Mhos/cm.                          |
| Resolution               | Mhos : 0.1 Micro mhos/cm  |
| Measuring Frequency      | 1000 HZ   |
| Temperature Compensation | 0 to 50°C   |
| Cell Constant            | 0.4 to 1.5 adjustable   |
| Function Selector        | CHECK/TDS./ COND./CELL CONST.                                       |
| CAL Function             | For full scale adj. (at back panel) Select at CHECK Position        |
| Digital Display          | 3 1/2 digit LED display   |
| Power Supply             | 220V $\pm$ 10% Hz AC  |
| Dimension                | L275 x B 175 x H 75 mm  |
| Weight                   | 2 Kg  |



**LABTRONICS**<sup>®</sup>  
Always a step ahead in technology...

All Designs And Specifications are subject to change without prior notice.