

# LT-27/ LT-18 DIGITAL D.O METER

## PRODUCT FEATURES

- ✓ D.O and Temperature indication (LT-27)
- ✓ Bright LED display
- ✓ Compact design
- ✓ Accurate and stable result

## ACCESSORIES

- ✓ Gold/Silver D.O Probe
- ✓ D.O. Membrane
- ✓ Temperature probe (LT-27)
- ✓ Rotational stand
- ✓ Instruction manual
- ✓ Dust cover
- ✓ Instruction manual



## TECHNICAL SPECIFICATIONS

Range	0 to 20.0 ppm
Temperature	0 to 100 (in model LT-27)
Resolution	0.1 ppm
Accuracy	0.1ppm $\pm$ 1 count
Temp. Compensation	0°C to 50°C manual
Compensation	$\pm$ 3%
D.O. Sensor	Amperometric gold/silver membrane type.
Display	3½ digit LED display
Power	230V A.C 50 HZ
Weight	2 Kg