

LT-32 DIGITAL POTENTIO- METER

PRODUCT FEATURES

- ✓ Highly stable & accurate
- ✓ 3½ digit LED display
- ✓ Rugged & economical
- ✓ Compact Design

ACCESSORIES

- ✓ Calomel electrode
- ✓ Glass electrode
- ✓ Platinum electrode (optional)
- ✓ Silver electrode (optional)
- ✓ Rotational electrode stand



TECHNICAL SPECIFICATIONS

Range	0 – 1999mV
Resolution	1 mV
Accuracy	± 1 mV
Input Impedancy	1012 Ohms.
Operating Temperature	10°C to 50°C
Display	3½ Digit 0.5 seven segment LED Indication with auto polarity indication
Power	230V ± 10 % AC, 50Hz
Dimensions	76x275x170 mm
Weight	2 kg. (approx)



LABTRONICS®
Always a step ahead in technology...

All Designs And Specifications are subject to change without prior notice.