

# LT-33/ LT-34 DIGITAL TURBIDITY METERS

## PRODUCT FEATURES

- ✓ Display - 3½ digit LED
- ✓ 30mm flat bottom Test Tube
- ✓ Range 0-1000 NTU (in Model 33)
- ✓ Ranges 0-200 NTU (in Model 34)
- ✓ Highly Stable & Accurate

## ACCESSORIES

- ✓ Test Tube 30mm dia 1 Set
- ✓ Instruction Manual 1 No.
- ✓ Dust Cover 1 No



## TECHNICAL SPECIFICATIONS

Product	Digital Turbidity Metter (LT-33)	Digital Nephelometer (LT-34)
Display	3½ Digit LED	3½ Digit LED
Detector	Photocell/Photodiode	Photocell/Photodiode
Range	0 to 1000 NTU in 2 ranges	0 to 200 NTU in 2 ranges
Resolution	1 NTU	0.1 NTU
ACCURACY	± 3% FS ± 1digit	± 2% FS ± 1digit
Sample System	60 x 25mm Clear Glass Test Tube	
Light Source	LED	
Power	230V ± 10%, 50Hz. A	
Dimensions	280 x 175 x 110 mm	
Weight	3kg approx	