

LT-331 DIGITAL NEPHELO TURBIDITY METER

PRODUCT FEATURES

- ✓ Display - 3½ digit LED
- ✓ 30mm clear glass Test Tube
- ✓ 3 ranges (0 to 20 NTU/JTU), (0 to 200 NTU/JTU), (200 to 1000 NTU/JTU)
- ✓ Highly Stable & Accurate

ACCESSORIES

- ✓ Test Tube 30mm dia 1 Set
- ✓ Instruction Manual 1 No.
- ✓ Dust Cover 1 No



TECHNICAL SPECIFICATIONS

Range	3 ranges (0 to 20 NTU/JTU), (0 to 200 NTU/JTU) (200 to 1000 NTU/JTU)
Accuracy	± 3% of full scale deflection in 0 to 200 NTU/JTU
Test Tube System	30 mm clear glass test tube.
Light Source:	6 V, 1 Amp. Tungsten lamp.
Display	3½ digit LED display
Detector	Photocell/ Photodiode
Calibration	With formazine solution
Power Supply	230 V ± 10% AC, 50Hz
Dimensions	(275L x 195W x 105H)mm
Weight	3kg