

LT-350 ULTRASONIC CLEANERS



PRODUCT FEATURES

- ✓ Stainless steel tank has resistance to wear and long work life.
- ✓ With digital timer & heater control.
- ✓ With Stainless steel basket
- ✓ Use just tap water, or industrial alcohol and solvent cleaner for higher cleaning requirement
- ✓ industrial control chip micro controller flexible circuit boards control
- ✓ Secure & stable
- ✓ Temperature tunable from 20°C to 80°C.

ACCESSORIES

- ✓ SS Basket
- ✓ Main Lead
- ✓ Instruction Manual



TECHNICAL SPECIFICATIONS

Ultrasonic Frequency	40,000 Hz
Material	Stainless Steel SUS304
Tank Capacity	2.0L, 3.5L, 4.5L, 6L
Timer	1 - 30 minutes, Digital timer
Power Supply	AC 100 - 120V, 50 / 60Hz AC 220 - 240V 50 / 60 Hz 240W 200 W,
Ultrasonic Power	240W
Heating Power	200 W, Digital Heater
Unit Size	325 x 265 x 280 mm
Tank Size	300 x 240 x 150 mm
Package Size	430 x 365 x 420 mm
Net. Weight	6.0 kg
Gross. Weight	7.0 kg