

LT-3510

DIGITAL PROBE SONICATOR

SOUND PROOF ENCLOSURE
(IN-BUILT)



PRODUCT FEATURES

- ✓ Microprocessor controlled & completely programmable
- ✓ Large touch screen display
- ✓ Temperature indicator & max temp. alarm
- ✓ Motorized operation table for ease
- ✓ Auto-frequency for convenience & optimal processing
- ✓ Variable power setting
- ✓ Auto frequency tracking
- ✓ Safety Alarms

ACCESSORIES

- ✓ 6 mm Probe
- ✓ Power cable
- ✓ Instruction Manual



ACCURATE
& STABLE



TOUCH
SCREEN



SOFTWARE
CONTROL SYS.



EASY
OPERATION



24
HOURS
WORKING

TECHNICAL SPECIFICATIONS

Frequency	18-28 kHz Auto Tracking
Frequency tracking accuracy	1 Hz real-time tracking
Processing Capacity	0.5ml to 500 ml
Standard Probe	Ø 6 mm (1/4")
Display	8" Touch screen
Display Function	Temperature, Power, Time, Gap Time
Temperature Range	0-99°C
Timer	Yes upto 99 hours, 59 mins, 59 secs
Alarm Functions	Over temperature, over current, over load, transducer abnormality, etc.
Operating Table	Electrically operated motorized
Lighting	Available inside cabinet
Data Memory	9 Memory Group
Output Power	750W (adjustable)
Dimensions	280 x 280 x 490 mm / 10.2 kgs



LABTRONICS®
Always a step ahead in technology...

All Designs And Specifications are subject to change without prior notice.