

# LT-37/ LT-36 MICROPROCESSOR COLONY COUNTERS

## PRODUCT FEATURES

- ✓ Count Correction Facility
- ✓ Result-Count Storage
- ✓ Audible Confirmation of Each Count
- ✓ Auto Marker Pen
- ✓ Uniform Glare Free Illumination
- ✓ Wolfhugel Glass Grid with Focussing Facility
- ✓ LED Based Illumination

## ACCESSORIES

- ✓ Counting Pen
- ✓ Magnifier Glass
- ✓ Stand
- ✓ Instruction Manual



LED DISPLAY



ACCURATE & STABLE



ACCURATE COUNTING



LED BASED ILLUMINATION



EASY OPERATION

## TECHNICAL SPECIFICATIONS

Display	3- Digit Seven Segment LED in Model LT-36 4- Digit Seven Segment LED in Model LT-37
Max Count	999 in Model LT-36, 9999 in Model LT-37
Keyboard	6 soft touch membrane keys
Correction	Count correction facility available to cancel the unwanted counts
Storage	40 Result-count storage retained even if the system is switched off
Average	Average mode facility available
Dish Size	110 mm
Magnification	x 1.7 surface magnification
Dimensions (I X B X H)	255 x 220 x 100 mm (Approx.)
Weight	3 Kg, (Approx.)
Power	230 V ± 10% AC, 50 Hz



**LABTRONICS®**  
Always a step ahead in technology...

All Designs And Specifications are subject to change without prior notice.