

LT-376 MICROPROCESSOR COLONY COUNTERS

6 DIGIT

PRODUCT FEATURES

- ✓ Count Correction Facility
- ✓ Result-Count Storage
- ✓ Audible Confirmation of Each Count
- ✓ Auto Marker Pen
- ✓ Uniform Glare Free Illumination
- ✓ Wolffhugel Glass Grid with Focussing Facility
- ✓ LED Based Illumination

ACCESSORIES

- ✓ Counting Pen
- ✓ Magnifier Glass
- ✓ Stand
- ✓ Instruction Manual



TECHNICAL SPECIFICATIONS

Display	6- Digit Seven Segment LED
Max Count	999999
Keyboard	6 soft touch membrane keys
Correction	Count correction facility available to cancel the unwanted counts
Average	Average mode facility available
Dish Size	110 mm
Magnification	x1.7 surface magnification
Reset	Data reset facility available
Memory	99 memory location
Storage	99 result count storage retained even if the system is switched off
Dimensions (L x B x H)	255 x 220 x 100 mm (Approx.)
Weight	3 Kg, (Approx.)
Power	230 V ± 10% AC, 50 Hz