

LT-42 DIGITAL DRY BATH INCUBATOR

PRODUCT FEATURES

- ✓ Uniform and accurate heating
- ✓ Over heat protection
- ✓ Available with buzzer alarm
- ✓ Highly stable and accurate
- ✓ Rock solid Block

ACCESSORIES

- ✓ Main lead
- ✓ Instruction manual



TECHNICAL SPECIFICATIONS

Display	16 X 2 Line Alpha Numeric LCD
Display Format	Menu driven
Setting	By Soft Touch Keys
Heating Devices	Test Tube Dry Bath
Temp. Sensing Devices	RTD/ PT100/ Digital Semiconductor
Range	5°C above ambient to 70°C
Resolution	1°C
Accuracy	± 1°C
Calibration	Manual (if required)
Block size	24 x 12 mm Tube holders.
Block Material	Rock solid Aluminum block.
Power Supply	230 V ±10% AC 50 Hz
Weight	3kg