

LT-422

DIGITAL DRY BATH INCUBATOR

PRODUCT FEATURES

- ✓ Uniform and accurate heating
- ✓ Over heat protection
- ✓ Available with buzzer alarm
- ✓ Highly stable and accurate
- ✓ Rock Solid Block
- ✓ Temperature Range Upto 150 °C

ACCESSORIES

- ✓ Main lead
- ✓ Instruction manual



TECHNICAL SPECIFICATIONS

Display	16 X 2 Line Alpha Numeric LCD
Display Format	Menu driven
Heating Block	24 Holes for 12.5 mm tubes
Setting	By Soft Touch Keys
Heating Devices	Test Tube Dry Bath
Temp. Sensing Devices	RTD/ PT100/ Digital Semiconductor
Range	5°C above ambient to 150°C
Resolution	1°C
Accuracy	± 1°C
Calibration	Manual (if required)
Block size	24 x 12 mm Tube holders.
Block Material	Rock solid Aluminum block.
Power Supply	230 V ± 10% AC 50 HZ
Weight	3kg