

LT-430 DIGITAL BLOOD ROLLER MIXER

PRODUCT FEATURES

- ✓ Rocking and rolling action for complete mixing
- ✓ six roller design with small footprint
- ✓ designed for continuous quiet operation
- ✓ adjustable speed range from 10 to 70 rpm

ACCESSORIES

- ✓ Main lead
- ✓ Dust cover
- ✓ Instruction manual



TECHNICAL SPECIFICATIONS

Display	LCD
Adjustable speed range	10 to 70 rpm
No. of Roller	6
Maximum load	10 test tube
Timer	Yes
Speed	Adjustable
Weight	3kg (approx)
Power supply	230V ± 10% ac 50Hz