

LT-44 ROTARY BLOOD MIXER



PRODUCT FEATURES

- ✓ AC Synchronous motor for noiseless operation.
- ✓ Suitable for a variety of centrifugal tube 1.5-50ml.
- ✓ Circular rotating action gives gentle and effective mixing.
- ✓ Fix speed of 30 RPM with an accuracy of ± 2 RPM.
- ✓ Continuous/infinite mode of operation.
- ✓ Fix plate angle of 45 Degree.
- ✓ Easy to operate, long running, stable and low noise.
- ✓ Maximum operational weight capacity of 1.5 kgs.
- ✓ Easy to remove and install disk plates for replacement and cleaning.

TECHNICAL SPECIFICATIONS

Speed	30 RPM
Slope angle	45°
Voltage	220 V AC
Power	6V
Maximum operational weight capacity	1.5 kg
Dimension(WDXH) MM	165X150X145



ACCURATE & STABLE



LOW NOISE



FIXED SPEED
30 RPM



EASY
OPERATION



24
HOURS
WORKING

LT- 43 BLOOD ROLLER MIXER



PRODUCT FEATURES

- ✓ Rocking and rolling action for complete Mixing
- ✓ Four roller design with small footprint
- ✓ Designed for continuous quiet operation

TECHNICAL SPECIFICATIONS

Roller Size	230mm
Controls	Analog
No. Of Rollers	4
Speed	33rpm (Approx.)
Maximum Load	5 test tubes
Timer	No
Power Supply	230V \pm 10% AC 50Hz



FIXED SPEED
25-30 RPM



ACCURATE & STABLE



EASY
OPERATION



LOW NOISE



24
HOURS
WORKING



LABTRONICS®
Always a step ahead in technology...

All Designs And Specifications are subject to change without prior notice.