

LT-47/ LT-46 MINI CENTRIFUGE

Mini Centrifuge equipped with a variety of centrifuge rotors and a variety of test tube sets. With flip open function, automatically stop when the cover is opened. It can be used for serum extraction from whole blood, various sample extraction, rapid sample sedimentation, separation of trace blood cells, microbial sample processing, PCR experiment zone centrifugation, etc.





PRODUCT FEATURES

- ✓ The unique rotor can meet different experiments without change
- ✓ Micro DC motor and carbon fiber rotor
- ✓ With silent technology operate with ultra-low noise
- ✓ Strong shell made of high-strength transparent material, resistance to fall and shock
- ✓ Flip switch function, automatically stop when the cover is opened



TECHNICAL SPECIFICATIONS

Image		
Product Model	LT-47	LT-46
Digital Display	Yes	No
Speed Range RPM	6000	6000
Controller Type	Adjustable speed & Timing	Fixed speed
Size	180x118x160mm	
Power Consumption	9w	
Input Voltage	24vdc 1.5 amp	
Motor type	Dc Permanent Magnet Motor	
Weight	500 gm	