

# LT-52 MICROPROCESSOR D.O METER

## PRODUCT FEATURES

- ✓ Data storage facilities
- ✓ Compact design
- ✓ Computer interface available
- ✓ Gold/Silver D.O probe
- ✓ Accurate and stable results

## ACCESSORIES

- ✓ Gold/silver D.O Probe
- ✓ D.O membrane
- ✓ Rotational D.O probe stand
- ✓ Temperature probe
- ✓ Dust cover
- ✓ Instruction manual



## TECHNICAL SPECIFICATIONS

Display	16x2 LCD display
Display format	Temperature and D.O simultaneously
Storage	1000 sample with date and time
Interface	USB
Do Range	0 To 20 ppm
Resolution	0.1 ppm
Temperature Comp	Automatic With Pt 100 Probe & Manual With Keys
Pressure Adjustment	0 To 999 Mm 01 Hg manual
Salinity Correction	0 To 20 Grill (manual)
Calibration	Automatic
Accuracy	± 2 % +1 Digit1
Temperature Range	-20.0 TO +200°C
Resolution	0.1°C
Accuracy	±0.2 %
Input Type	Pt 100
Calibration	Manual
Power	230V A.C 50 HZ
Weight	1 Kg