

LT-63 AUTO KARL FISCHER TITRIMETER

PRODUCT FEATURES

- ✓ Auto Titration Control
- ✓ Auto Zero Burette
- ✓ Adjustable Timer
- ✓ Visual & Audible Alarm

ACCESSORIES

- ✓ Dual platinum electrode
- ✓ Burette clear glass 10ml
- ✓ Rubber bellow
- ✓ Titration Vessel
- ✓ Reagent Bottel



LED DISPLAY



INBUILT
STIRRER



STABLE
DESIGN



ACCURATE
& STABLE



EASY
OPERATION

TECHNICAL SPECIFICATIONS

Display	Array of 10 LEDS
Sample Size	1-50mg of H ₂ O
Indication	RUN, ADD, END LEDS
Kf Flow	Auto Control Valve
Burette	Auto Zero 10ml Burette
Timer	Adjustable Timer
Dispensing Resolution	0.05ml
Measuring Electrode	Dual Platinum KF Electrode
Power	230V ± 10%, 50Hz. AC