

LT-70/ LT-71 DIGITAL FRIABILITY TEST APPARATUS

PRODUCT FEATURES

- ✓ Programmable Time With Digital Display
- ✓ Auto Stop After Completion of Set Time With Alarm
- ✓ Auto Removal of Sample Without Opening the Drum
- ✓ Precise Control of RPM 25 ± 1

ACCESSORIES

- ✓ Acrylic Drum, 2Nos (1Nos acrylic drum in model LT-71)
- ✓ Main lead
- ✓ Instruction manual



TECHNICAL SPECIFICATIONS

		
	LT-71	LT-70
Model		
Rotation Speed	25 ± 1	
Time Range	1 sec to 99 min 59 sec	
No Of Wheels	ONE	TWO
Test Time Programmable	Programmable	
Electric Motor	Low rpm A.C. Motor	
Auto Removal of sample	Yes	
Audio Alarm at completion	Yes	
Power	230V ± 10% AC 55hz	
Dimension	270x280x290mm (approx)	
Weight	10kg	