

# LT-730 DIGITAL AUTOMATIC VISCOSITY METER

## PRODUCT FEATURES

- ✓ Stepping Motor for Accurate, reliable operation
- ✓ Direct readout of all measurement parameters
- ✓ Auto Ranging
- ✓ Time Function for measurement
- ✓ Accurate calibration

## ACCESSORIES

- ✓ Spindle 4Nos
- ✓ Instrument stand 1Nos
- ✓ RS-232 cable 1Nos
- ✓ Instruction manual
- ✓ Power adapter



LCD  
DISPLAY



STABLE  
DESIGN



ACCURATE  
& STABLE



EASY  
OPERATION



RS232  
CONNECTIVITY

## TECHNICAL SPECIFICATIONS

Measuring range	10 – 2,000,000 C.P
Display Resolution	0.01 C.P
Measurement accuracy	±5% Newtonian Fluid)
Speed	0.3 0.6 1.5 3 6 12 30 60 RPM/min
Outline size	95*130*155 Excluding base Net weight 2kg Excluding base
No of standard spindel	#1, #2, #3, #4 is standard (#0 is optional )
Display	LCD display with backlit
Power Supply	230VAC ±10 %, 50Hz
Dimension	400x350x360mm approx
Weight	9kg approx