

LT-7310 TOUCH SCREEN VISCOMETER

PRODUCT FEATURES

- ✓ Stepping Motor for Accurate, reliable operation
- ✓ Direct readout of all measurement parameters
- ✓ Auto Ranging
- ✓ Time Function for measurement
- ✓ Sound alarm at under 20% Torque
- ✓ Accurate calibration

ACCESSORIES

- ✓ Spindle 4Nos
- ✓ Instrument stand 1Nos
- ✓ Touch Pen
- ✓ RS-232 cable 1Nos
- ✓ Instruction manual
- ✓ Power adapter



TECHNICAL SPECIFICATIONS

Measuring Range (mPa.s)	20-2,000,000 (mPa.s)
RPM	0.3, 0.6, 1.5, 3, 6, 12, 30, 60
No of Standard Spindles	#1, #2, #3, #4 is standard (#0 is optional)
Measurement Accuracy	± 1.0% of range
Repeatability	± 0.5%
Display	LCD Touch Screen
Temperature RTD Monitoring Sensor Probe	Yes
Output	RS 232 Interface
Power Supply	AC 220V/50Hz
Dimension LxWxH (mm)	400 x 350 x 360
Weight	9 Kg