

# LT-74/ LT-73 DIGITAL DISINTEGRATION TEST APPARATUS

## PRODUCT FEATURES

- ✓ Programmable Time With Digital Display
- ✓ Auto Stop After Completion of Set Time With Alarm
- ✓ Memory for Protection of Running Test in Case of Power Failure
- ✓ Separate Regulator for Control of Heating
- ✓ Facility for Hold and Resume the Running Test for inspection of Sample

## ACCESSORIES

- ✓ External Temperature Sensor
- ✓ Mains Lead
- ✓ Instruction Manual
- ✓ 1000 ml Glass Beaker 2 Nos.
- ✓ Disc 12 Nos.
- ✓ Tube 12 Nos.



## TECHNICAL SPECIFICATIONS

Type of Heating	Hot Plate
Up-down Speed	30 ± 2 strokes/min
Up-down Length	55 ± 5mm
Time Control	Digital (Programmable from 1 sec to 99hr, 59min, 59 sec)
Electric Motor	Low rpm A.C. Motor
No. of Basket	Two (One in Model 73)
Audio Alarm at completion	Yes
Hold & Resume Facility	Yes
Reaction beaker	1000ml
Resolution	±2%
Power	230V ±10% AC 50hz
Dimension	365x350x495mm (approx)
Weight	9kg (approx)