

LT-77 BULK DENSITY APPARATUS

PRODUCT FEATURES

- ✓ Extremely useful for Pharmaceutical Lab
- ✓ Highly Accurate and Easy to Operate
- ✓ Microcontroller Based Digital Display
- ✓ Soft Touch Membrane Keys
- ✓ 3 LED's for status Indication
- ✓ Audible Alarm Sound
- ✓ Last Count Retention on Power Failure

ACCESSORIES

- ✓ Glass cylinder 2 Nos(100ml or 250ml)
- ✓ Cylinder holder
- ✓ Main lead
- ✓ Instruction manual



TECHNICAL SPECIFICATIONS

Maximum Counts	9999
Display	4 Digit 7-Segment Red LED
Keyboard	6 Soft Touch Membrane Keys
Cylinder	2 Nos., 100ml/250 ml Cylinders
Capacity of cylinder	100ml/250ml(mention during orders)
Drop Height	25 mm (Approx.)
Power	230 V \pm 10% AC, 50 Hz.
Dimensions	600 x 450 x 200 mm (Approx.)
Weight	12 Kg (approx.)