

LT-770 TOUCH SCREEN DIGITAL MOISTURE BALANCE

IR SOURCE

PRODUCT FEATURES

- ✓ Long life IR based heating
- ✓ Highly stable and accurate
- ✓ Colored touch screen
- ✓ Fast heating

ACCESSORIES

- ✓ Touch Pen
- ✓ Sample Tray
- ✓ Main lead
- ✓ Instruction manual
- ✓ Weight box



IR Source



TOUCH SCREEN



SOFTWARE CONTROL SYS.



ACCURATE & STABLE



EASY OPERATION



FAST HEATING

TECHNICAL SPECIFICATIONS

Max. Capacity	120g
Accuracy	0.005g
Readability	0.01%
Weighing sensor	HBM
Temperature range	40-230°
Heating source	High efficiency IR source
Heating way	Gentle, standard, fast
Halt mode	Automatic, manual, timed
Result display	Test temperature, test time, moisture content, solid content, dry weight, dynamic curve, moisture regain, solid content, weight loss, dry weight, test method, sample name
Tray material	Stainless
Tray diameter	Φ90mm
Interface	Standard RS232 communication interface, which can be connected to printers, PCs and other peripheral devices, and USB can be customized
User management	Password protection, multiple users, auto-lock
Display Screen	5"HD touch screen
Product size	370x215x195mm
Package size	445x300x330mm
Weight	5kg approx



LABTRONICS®
Always a step ahead in technology...

All Designs And Specifications are subject to change without prior notice.