

LT-780 DIGITAL TABLET HARDNESS TESTER




TECHNICAL SPECIFICATION

Display	LCD with backlit display
Keyboard	User Friendly 9 Keys operations
Hardness Range	5-500 N
Hardness Accuracy	± 0.05% FS
Measuring Units	N(Newton), KGF (Kilogram Force) and Lbf (Low bound Force)
Interface Ports	RS-232
Dimensions	360 x 65 x 100 mm
Power	12V , 300 mA DC power adapter
Battery	1200 mAH, 8.4V, Ni-Hi Rechargeable Battery
Peak Force Measurement	Yes
Acrylic Safety Glass Cover	Yes
Stored Results	Upto 10 (also calculate the average value of test results)
Tablet Size	Yes (2mm to 28 mm)



LT-113 DIGITAL HAEMOGLOBIN METER



ACCESSORIES

- ✓ Glass test tubes 5 Nos.
- ✓ Instruments cover 1 Nos.
- ✓ Operating Manual 1 Nos.
- ✓ Spare Lamp 1 Nos. (Fitted inside)

TECHNICAL SPECIFICATIONS

Display	3 Digit Red LED Display.
Display Modes	% T, O.D. and 3 Concentration Modes C1,C2,C3
Zero Adjustment	Automatic (soft touch key)
Accuracy	± 0.02 O.D.
Stability	± 0.02 O.D. Per Hour.
Light Source	3.5V LED bulb
Amplifier	Micro Controller Based.
Memory	In built memory to retain last "K" factor for concentration even after power loss.
Sample Volume	1 ml. Minimum.
Power Supply	230 V ± 10 %, 50 Hz.