

# LT-84/ LT-83

## MOTORLESS MAGNETIC STIRRERS



### PRODUCT FEATURES

- ✓ Digital display provided easy reading of the speed and time.
- ✓ The most advance magnetic coil technology for wear-free drive without moving parts.
- ✓ Automatic reverse rotation (Pulse function) every 30 second provided excellent mixing performance.
- ✓ Max stirring quantity(H<sub>2</sub>O) of 800 ml.
- ✓ Speed range of 15 to 1500 rpm.
- ✓ Slim design with unique flat surface design.
- ✓ Water spill resistant design – IP 65.



### ACCESSORIES

- ✓ Model: LT-84. Stirrer bar 2 piece, 9V 1.5 power adaptor, operation manual
- ✓ Model: LT-83. Stirrer bar 1 piece, 9V 1.5 power adaptor, operation manual

### TECHNICAL SPECIFICATIONS

Image		
Product Model	LT-84	LT-83
Speed Limit	15-1500 RPM	
Stirring Position	2	1
increment and decrements	50 RPM	
Display	Digital LED	
Working Plate Material	SS Plate	
Setting Dia.	105 mm	
Maximum Stirring Capacity	1000 ml H <sub>2</sub> O	
Recommended Stirrer Bar Size	25x8mm PTFE Coated	
Pulse Mode	30 second auto reverse direction	
Timer	99 Minuits and Inifinité Mode	
Last run Memory Save	Yes	
Protection Class	Ip65	
Permissible Ambient Temp	5 °C to 40 °C	
Permissible relative Humidity	80%	
Voltage	100-240 VAC	
Frequency	50/60 Hz	
Weight	0.5 Kg	0.3 Kg
Dimension(WxDxH)	230x12x175 mm	162x12x110 mm