

LT-91/90/89

MULTI POSITION

MAGNETIC

STIRRERS






Multi position magnetic stiner adopt international advanced design concept and manufacturing technology, with smooth and beautiful appearance, compact and firm structure, high universality, simple operation and other advantages. It is a rare instrument for scientific research and experiment. This magnetic stirrer is mainly applicable to laboratory instruments for liquid mixing, mainly used for stirring or heating and stirring low viscosity liquid or solid-liquid mixture. It is widely used as an experimental tool in colleges, environmental protection, sanitation, epidemic prevention, chemical.

PRODUCT FEATURES

- ✓ High synchronization rate, speed error of each point < 10rpm
- ✓ LCD display, the display content is beautiful and clear, and the reading is accurate and intuitive
- ✓ Specially made film coated disk surface, which is anti-skid, scratch resistant, corrosion resistant and easy to locate the motor
- ✓ Push type knob, long press to start operation, sensitive response
- ✓ The Rubber foot has a stable structure to prevent the instrument from sliding and toppling.



TECHNICAL SPECIFICATIONS

Image			
Product Model	LT-91	LT-90	LT-89
Voltage V	100-240		
Frequency Hz	50-60		
Power W	135	90	45
Display Method	LCD		
Timing Range	1-999min		
Maximum Stirring Volume(H ₂ O) ml	1000		
Stirring Points	9	6	3
Speed Range rpm	100-1500		
Speed Control Accuracy rpm	±10		
Motor Type	DC Brushless Motor		
Working Plate Size mm	390*390	260*390	190*390