






LT-14/17/19/28

PORTABLE INSTRUMENTS



TECHNICAL SPECIFICATIONS

Model No.				
Model No.	LT-14	LT-17	LT-19	LT-28
Range	0-14 pH	0-1999 μ S, 0.19.99 μ S	0.1999ppm, 0-19.99 ppt	0-19.9 ppm
Resolution	0.01 pH	1 μ S	1 ppm	0.1 ppm
Impedance	>10 ¹² ohms	-	-	-
Sensor	combination pH electrode	Cond. Cell (dip type)	TDS cell (Dip type)	Amperometric D.O probe
Battery	9V DC	9V DC	9V DC	9V DC
Accuracy	0.01 pH	\pm 1% FS	\pm 1% FS	\pm 2% FS
Temp. Comp.	0-100° C (Manual)	-	-	5-55° C
Display	3½ digit LCD	3½ digit LCD	3½ digit LCD	3½ digit LCD
Size	180 x 85x 30 mm	180 x 85x 30 mm	180 x 85x 30 mm	180 x 85x 30 mm
Weight	250 gm (Approx)	250 gm (Approx)	250 gm (Approx)	250 gm.(Approx)